

# 3 Compartment Sink

## SINK 1 – WASH

- Warm water and detergent

- Sort, scrape, and pre-rinse

## SINK 2 – RINSE

- Clean water at a minimum temperature of 43°C (110°F)

## SINK 3 – SANITIZE

- Sanitize for at least 45 seconds in any one of these:
  - Chlorine\* at 100 ppm
  - Quaternary ammonia at 200 ppm
  - Iodine at 25 ppm
  - Acceptable sanitizer
  - Clean water at a minimum temperature of 77°C (170°F)

\*Use 2 mL of chlorine (bleach) for every 1L of water to yield ~ 100 ppm solution

- Air dry

